IT

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## INDICATIONS FOR USE

The prosthetic components are designed and manufactured to allow correct and optimal screwing of the prosthesis onto the various implant components already on the market.

### V4 BASE - CEMENT-RETAINED TECHNIQUE CONVERTERS

The cement-retained technique converters are made of grade 5 titanium, both with non-rotational engagement (for single prosthesis) and with rotational engagement (for multiple prostheses, Implant Bridge Toronto); this type of product can be modified by the dental technician only in the height of the cementable cannula, as required, through an accessory made specifically to optimise this operation. The cement-retained technique converter is sold with a single permanent prosthetic screw produced by New Ancorvis. The types of converters are produced to be used as both temporary and permanent prostheses.

#### Technical specifications:

- Grade 5 titanium to ensure excellent resistance to the biting and chewing stress.
- Available with transmucosal heights to better adapt to the patient's aesthetic requirements. Observe the reference table.
- To connect to the implant, use the screw included in the package.

### VITI HIGH GRI

The V4 Base is sold with a permanent prosthetic screw. It should only be used to fix the prosthesis in the oral cavity; the purpose of the screws is to permanently tighten the prosthesis to the implant; they are characterised by a special surface coating to reduce the unscrewing phenomenon. The prosthetic screws made by New Ancorvis are compatible and therefore replaceable with the original screws of the implant manufacturer.

#### Technical specifications

- The screws are made of grade 5 titanium for excellent resistance to biting and chewing stress.
- The screws are included inside the converter package.
- The coating applied to the screws is: Balinit C 1500 (composition: a-C: H: W, colour: anthracite).

### V4 LEGO - SUPPORT FOR REINFORCED PROVISIONAL PROSTHESIS

V4 Lego is a support for reinforced provisional prostheses to be completed with a surface polymer coating. It can be made in either Gr.5 Titanium or cobalt chromium using the selective laser melting technique. V4 Lego is interfaced with the V4 base thanks to a special threading that allows height adjustment and orientation in the space, based on the planned prosthetic reconstruction. Once the height has been adjusted, it can be secured through the cement-retained technique using the dedicated hole to inject a retaining cement (for the choice of material see the indications provided by the manufacturer of the cement). To bind the various components, it is advisable to use a medical grade steel wire/strap/bar. Any changes made by the user to the structure of the MD may compromise its correct operation.

#### Technical specifications:

- Available in both grade 5 titanium and cobalt chromium to guarantee excellent resistance to the biting and chewing stress produced in SLM.
- Available in different shapes to adapt to various situations.

### V4 PICK UP

V4 Pick Up is a special transfer that can only be used with the V4 base line; these components are distinguished by laser screen printing and by anodisation to identify the family they belong to (according to the implant house) in order to correctly couple the V4 base and the V4 Pick Up. It is advisable to check the correct tightening between V4 Base and V4 Pick Up in order to guarantee a correct impression (the use of an enlargement system is recommended). Ensure the transfer is tightened until contact with the cement-retained converter). For use in the impression coping configuration, the dedicated coping screw must be used for each implant platform; all coping screws have a knurled part to facilitate manual tightening and engagement for a 1.27mm hex driver (New Ancorvis coding: B).

#### Technical specifications

- The V4 Pick Up is made of grade 5 titanium.
- V4 Coping screw produced in grade 5 titanium; for tightening use 1.27mm hex driver (New Ancorvis code: B).
- Laser screen printing and anodisation to identify the various implant families. See the reference table.

### V4 SCAN MARKER

The V4 Scan Marker is the scan-abutment made exclusively for the V4 Base, therefore use is linked only to this type of product. The insertion between the two components is unique in that the Scan Marker can be inserted in one position only, thanks to the three-lobed body of the cement-retained converter, ensuring a precise prosthetic restoration; the insertion phase is made easier as two notches have been made (one on the converter, one on the scan marker) which must be aligned. It is strictly necessary to use this scanning product only with the new V4 base component as, should it be necessary to adjust the height of the same, it cannot be digitally scanned. Vertical positioning is guaranteed by a click, which certifies correct positioning in "Z"; if after several uses the "click" can no longer be heard, the V4 Scan Marker should be replaced. It is advisable to keep the part inclined buccally.

#### Technical specifications

- The V4 Scan Marker is made of Peek.
- Product to be used only on the V4 Base.
- Use the specific libraries for this product.
- In the case of creation of Implant Bridge, use this product for no more than two implant elements; for larger workings use the AQ scan abutment.

## WARNINGS

- The screws, the screw housings and the platforms must NEVER be tampered with.
- For screwing instructions, refer to the indications on the label.
- The product is disposable: any reuse exposes the user to the risk of cross-infection, therefore the appropriate reference symbols are affixed to the label and packaging.
- Store the product in a dry place at room temperature, taking care to keep the package intact.
- Dispose of it in accordance with current regulations.
- The dental technician is responsible for and takes on all the risks concerning:
  - The choice of the implant component compatible with the uses required by his/her activity;
  - The design of the overlying pillar using safe parameters and thicknesses to ensure precision and stability over time, relieving New Ancorvis srl of any liability linked to stability as it is an operation managed entirely by the dental technician in his/her laboratory.
- New Ancorvis shall not be held responsible for an direct or indirect damage caused by the improper use of the devices or tampering.
- All implant components must be used only by personnel with the necessary operating qualifications.
- We recommend reading the instructions for use of the implant to which the structure must be connected.
- The implant components are sold in closed plastic boxes with a label showing the manufacturer's data, the batch number, tightening torque where necessary, CE marking and harmonised symbols.

## CLEANSING, DISINFECTION AND STERILISATION

Correct cleaning, disinfection and sterilisation are essential to reduce and prevent risks.

To correct perform these operations, these instructions should be followed:

- The devices are supplied cleaned and washed according to the New Ancorvis protocol, suitably validated.
- The devices are supplied NON-STERILE, therefore before use on the patient it is mandatory to sterilise the device in a steam autoclave at 134°C for at least 5 minutes, according to the validated cycle). Before sterilisation, it is important that the device undergoes a cleaning and washing process with a method validated on the basis of the residues of the processes performed by the user.
- Higher temperatures may change the physical-chemical composition of the components. Do not use dry or chemical autoclaves.

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